

BINOCULAR ACCOMMODATION
AUTO REF/ KERATO METER
WAM-5500

ADVANCED ACCOMMODATION MEASUREMENT

Binocular, Dynamic, Refractive, Pupil Diameter Measurements

NEW

Specification

Refraction Measurement	Sphere -22~+22D (0.01/0.12/0.25DStep)
	Cylinder 0~±10D (0.01/0.12/0.25DStep)
	Axis angle 0~180° (1°Step)
Measurement of Corneal Radius	Corneal radius 5.0~10.0mm (0.01mmStep)
	Refractive power 33.75~67.5D (0.01/0.12/0.25DStep)
	Cylindrical power 0~±9D
	Axis angle 0~180° (1°Step)
Pupil Diameter Measurement	φ2~φ8 0.1mm step
Vertex Distance	0, 10, 12, 13.5, 15mm
Minimum Pupil Diameter	φ2.3mm
Pupillary Distance	Measurement Range 85mm (1mmStep)
Printer	Thermal printer with Automatic Cutter (Width 57mm)
Internal Monitor	5.6 inch LCD Display (Color)
Movable Distance	Back/Force ±17mm Right/Left ±43mm Up/Down ±15mm
Movable Distance of Chinrest	±30mm
Overall Dimension	(W) 327mm × (D) 496mm × (H) 515mm
Weight	About 20kg
Output	RS232C Interface
	Video Output (NTSC)
Rated Voltage	100~240V
	50/60Hz
Consumption	80VA
Power Save	OFF, 3, 5, 10 min. (Selectable)

- Specifications and design are subject to change without prior notice for improvement.
- The part of screen is composed photograph.
- The color of instrument on the catalogue and the real product might be different.
- Medical equipment manufacturer license No.37B3X00002


Accessories

- Test eye / 1 pc
- Power cable / 1 pc
- Printer paper / 3 rolls
- Fuse 2A / 2 pcs
- Chinrest paper / 1 pack (1000 sheets)
- Chinrest pin / 2 pcs
- Dust cover / 1 pc
- Fixation target / 1 pc
- Fogging lens / 1 pc
- Occluder / 1 pc
- Near point measurement unit / 1 set
- Contact lens holder / 1 pc

Option

- WCS-1 HI SPEED MODE Control software
- MDC-1 The Measurement Data Collection software
- Electric table
- Chair

Overall Dimension



Hardware Requirement
 OS : Japanese/English Version of Microsoft Windows95/98/2000/XP
 CPU : Pentium2 400Mhz or better (recommended)
 Memory : 128MB or more (recommended)
 (Mentioned company and product names are the trademark or the registered trademark of each company.)

HIGH-TECH
HIGH-TOUCH
HIGH-IDEALS
RyuSyo
 RyuSyo Industrial Co.,Ltd.
 958 Ikeuchi, Konan-cho, Kagawa-gun, Kagawa-ken 761-1494 Japan



CAUTION For the proper and safety use, please read the Instruction Manual thoroughly before using the product.

Grand Seiko Co.,Ltd.

Head Office & Plant
 12-19 Kashima, Fukuyama, Hiroshima 720-0091 Japan
 Tel.*81-84-952-2151 Fax.*81-84-951-9456
<http://www.grandseiko.com>

חברת מייצגת:

ד.ט.י. שחתי אופטיקה בע"מ
ציוד אופטלמי ורפואי




Grand Seiko Co.,Ltd.

ADVANCED ACCOMMODATION MEASUREMENT

Binocular, Dynamic, Refractive, Pupil Diameter Measurements

In addition to the Binocular Refractor and Keratometer measurement functions, the WAM-5500 can measure the Refraction and the Pupil Diameter simultaneously. These are advanced features provided on the earlier WR-5100K AutoRef/Keratometer version.

With this function, it can be used to confirm the change on the accommodation function before and after the surgery. In addition, it can be used to confirm presbyopia and the accommodation error due to fatigue as a result of VDT work or other eye fatigue instances. By connecting with the computer, you can measure refractive value at a rapid, continuous speed, thus allowing you to measure the accommodation response dynamically.



Selection of A.M.MODE


No.	00001
NAME	2005 3 9 14:30
VD=0	
<R>	FAR -4.75 -0.25 62
	-4.75 0.00 φ 5.0
	-4.75 0.00 φ 5.0
	-4.75 0.00 φ 5.0
	-4.75 0.00 φ 5.0
	-4.75 0.00 φ 5.0
<L>	-2.0D/ 50cm -4.50 -0.75 90
	-4.50 -0.75 φ 5.1
	-4.50 -0.75 φ 5.1
	-2.5D/ 40cm -4.50 -0.75 89
	-4.50 -0.75 φ 5.2
	-3.0 D/ 33cm -4.50 -0.75 90
	-4.50 -0.75 φ 5.1
PD = 64	
Grand Seiko Co., Ltd.	
WAM-5500 S/N:000001	

Right Eye Measurement Data

Left Eye Measurement Data

A.M. MODE

It measures the pupil diameter at 0.1mm step (minimum 2mm) simultaneously with refractive value. The measurements are then presented on the monitor. Subsequently, you can print the results and check the constriction of pupil. Studies have indicated that pupil constriction is evident when the eye is accommodating. Thus, the pupil diameter measurement function allows you to more precisely ascertain whether the subject's eyes are accommodating or not. Another option is to present the image via video thus giving you precise image analysis of the pupil constriction as well. The ideal refractive value for 50, 40, 33, 25, and 20cm is displayed, thus you can compare the patient data with the ideal value more easily as it relates to the actual distance of the Near Point Card.



Measurement with Both Eyes

With both eyes open, the natural condition, the measurement is taken. Thus, there would be less Instrumental Myopia, and allows you to achieve more precise measurement.

Easy Measurement with Open View Window


The Open View Window, which adopted the Half Mirror, makes it possible to take a measurement without the subject being aware as there is no distraction of view. For the infant/children, it gives less pressure, and you can measure easily and quickly. Additionally, it allows you to go through the measurement process smoothly since you can observe the condition of a subject directly.

Any Target

You can select from any target, as you like. Needless to say, it can be an eye-test chart. Or for the infant/children, you can choose from what they are interested in, such as toy and cartoon. It helps the measurement process go more smoothly. Additionally, since you can put a target at any distance, you can measure middle and near target as well.

HI SPEED MODE

By connecting with the computer and using our dedicated software - WCS-1 (optional), you can measure the SE value and Pupil Diameter at 0.2 seconds per step. In addition, you can analyze this data with spreadsheet software like Microsoft Excel. Unlike the previous model, WR-5100K, this Hi Speed Mode allows you to more dynamically measure the eye. Additionally, in this Mode, since the refractive value is calculated on the entire meridian instead of a specific one, you can understand the refractive value precisely.



● Example ; 33 years old male
This is the data of moving the Near Point Card every 20 seconds. With this data, the subject eyes are accommodating from the ideal value of -2.00D to -5.00D. Additionally, the pupil diameter is constricting from 4.0mm to 3.2mm. (Blinking at the same interval.)

Indication of SE (Spherical Equivalent)

You can confirm the change of accommodation as the SE is indicated on the monitor each time (3 data).

Easy Measurement

Since it is adopting the Open View Window, you can measure a subject with glasses and contact lens as well as a subject with the IOL, easily. It also allows you to confirm the prescription.



Data Output

Video Output / It allows to show the same image as internal monitor on an external monitor (NTSC System)
RS232C Interface / It allows you to send the measurement data.
Software / WCS-1 (Option) HI SPEED MODE Control software
MDC-1 (Option) The Measurement Data Collection software

Mode Change

K/R Mode Measures Ref and Kerato Value.
R, K Mode Measures Ref Value when "R" and Kerato Value when "K"
IOL Mode It makes easier to measure the IOL implanted subject. "IOL" is shown on the print out.
PD Measures the PD for far vision automatically, and calculate the PD for near vision (select the work distance from 30, 35, 40, 45cm).

Right eye data				Date & time	
No.	00001				
NAME	2005 3 9			14:30	
VD=12					
<R>	SPH	CYL	AX	Refractive data	Representative value
	-4.75	0.00			
	-4.75	0.00			
	-4.75	0.00			
	-4.75	0.00			
<R>	mm	D	AX	Corneal data	Representative value
R1	7.62	44.25			
R2	7.62	44.25			
AVE	7.62	44.25			
CYL		0.00			
<L>	SPH	CYL	AX	Refractive data	Representative value
	-4.75	0.00			
	-4.75	0.00			
	-4.75	0.00			
	-4.75	0.00			
<L>	mm	D	AX	Corneal data	Representative value
R1	7.62	44.25			
R2	7.62	44.25			
AVE	7.62	44.25			
CYL		0.00			
PD = 67 NPD=62					
Grand Seiko Co., Ltd.					
WAM-5500 S/N:000001					
Left eye data				Pupillary distance	

Left eye data

Pupillary distance